WIDE TACTICAL TOURNIQUET



Manufacturer code : 1401008 – Distributor code : S524284

Description : Tourniquet with a 4cm wide strap

The 4cm wide strap provides higher compression than most other standard turnstiles, which often have a wide buckle, but only 2.5cm strap; the comfort of the injured person is also improved.
It does not require any Velcro to maintain pressure, so it is more reliable. It is light and compact. The strap is strong. A label is sewn on the strap marked “TIME” to indicate the installation time.

Material : textile strap – fixing ring and hook steel-made
Coulor : black
Open dimensions : 3.8 X 105 cm
Folded dimensions: 4.5 X 15.5 cm
Weight : 0,107 Kg

Standards : Medical device class I according to Regulation (EU) 2017/745 CE

Instruction of use :
- Locate the tourniquet loop and pass it around the limb affected by the hemorrhage.
- Apply the tourniquet to the affected limb (arm or thigh):
  - above the elbow or knee,
  - between the wound and the heart,
  - never on a joint.
- Pull on the strap to tighten the buckle around the limb.
- Turn the aluminum handle to stop the bleeding.
- Lock the handle by inserting one end into the rigid plastic hook.
- Note the time of tourniquet application and monitor for bleeding to stop.
- Seek medical help immediately.

Precautions of use and contraindications
- For single use.
- Intended for external use.
- Keep out of reach of children.
- Do not use if sachet is opened or damaged.
- Store in a dry place and at room temperature.
- After use, throw away with medical waste.

Validity date : 5 years

Packing :

	Dimensions cm	Weight Kg	EAN code	Customs code
Individual pack	4,5x15,5 cm	0,107	6951454703860	3005109000
Carton of 100	44x28x32 cm	8,8		